

Содержание

3. Настройка сервиса для обогащения справочника GUI СКАТ данными о базовых станциях из внешнего источника по API	3
---	----------

3. Настройка сервиса для обогащения справочника GUI СКАТ данными о базовых станциях из внешнего источника по API

В случае отсутствия данных координат в справочнике для некоторых базовых станций предлагается использование **сервиса для обогащения справочника GUI СКАТ данными о базовых станциях из внешнего источника по API**.


Активация сервиса происходит добавлением в

```
/var/questor/backend/.env
```

Следующий параметров:

```
#Enables autoload of OpenCellID dictionaries
CELL_TOWER_AUTOLOAD_ENABLED=1
# 0 - Disabled (default)
# 1 - Enabled
#Token in OpenCellID platform
CELL_TOWER_LOADER_API_TOKEN=pk.fce9af82655c117ef42c88d36885cc8b
#Level of log
CELL_TOWER_LOG=1
# 0 - Info
# 1 - Debug
```

Получение токена осуществляется через портал <https://my.unwiredlabs.com/trial>. После регистрации токен будет доступен:



The screenshot shows the Unwired Labs Location API dashboard. The header includes the Unwired Labs logo, navigation links (Products, Use Cases, Pricing, Contact Us), and user links (Dashboard, Logout). The main content area is titled "Hi Vladislav!" and includes a welcome message. A sidebar on the left lists navigation options: API Sandbox, Reports, Devices, Geolocation API, Geocoding APIs, Maps, Account, Account details, API Access Tokens, and Logout. The main content area is titled "Manage your API Access Tokens" and includes a warning message: "Your plan allows 1 Access Token and you've reached this limit. To create more tokens, please upgrade your account." Below this is a table with one row showing an "Access Token 1" with a "Show Token" link and a "Created On" date of "26 August 2022 08:06 AM UTC".

Label	Access Token	Created On	
Access Token 1	Show Token	26 August 2022 08:06 AM UTC	View Logs

При нажатии *Show token* появится бесплатный ключ с ограничением 100 запросов в день.